ELANEE®

## Beckenboden-Trainingshilfen Phase I aktiv

<b>D</b>	Beckenboden-Trainingshilfe Phase I aktiv – Gebrauchsanweisung.....	2
<b>GB</b>	Pelvic Floor-Training Aids Phase I active – Instructions for Use.....	9
<b>F</b>	Accessoires de rééducation du plancher pelvien – Niveau I actif– Mode d'emploi.....	15
<b>I</b>	Ausilio per l'allenamento dei muscoli del diaframma pelvico fase I attiva – Istruzioni per l'uso.....	23
<b>NL</b>	Bekkenbodem-trainingshulpjes Fase I actief – Gebruiksaanwijzing.....	30
<b>TR</b>	Pelvik tabanı egzersiz yardımcı aşama I etkin – Kullanma kılavuzu.....	38
<b>RU</b>	конусы для укрепления мышц тазового дна Фаза I активная – Инструкция по применен.....	44

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### Instructions for use pelvic base training aids phase I

#### Dear Customer,

thank you very much for choosing our Pelvic Floor Training Aids Phase I active. These instructions for use will help you operate the product in the best possible way. Please read carefully through the instructions for use – in particular the “safety instructions” – prior to first usage.

Please contact us with any questions and suggestions!

*Your ELANEE-Team*

 Please observe both instructions for use and packaging and store them in a safe place for future reference!

*Made in Germany*

#### 1. Scope

These instructions for use apply to the following product:

**REF** 701-00 ELANEE Pelvic Floor Training Aids Phase I active

#### 2. Product content

ELANEE Pelvic Floor Training Aids Phase I active consists of a set of four tampon-shaped cones with a different colour and weight in each case:

<b>Colour</b>	yellow	blue	green	pink
<b>Weight</b>	20g	28g	52g	71g

Instructions for use are attached.

### 3. Purpose and intended use

Training with the Pelvic Floor Training Aids Phase I active ensures the prevention, reduction or elimination of stress incontinence caused by weakened pelvic floor muscles. The training is suitable

- during pregnancy and after birth for the prevention of weaknesses of the pelvic floor and to support regeneration and/or
- the reduction, elimination and prevention of urinary incontinence (as a result of load/stress incontinence) in women (especially during the climacteric period and menopause) .

The pelvic floor Training Aids Phase I active are not suitable for application purposes beyond the ones named.

### 4. Areas of application

- **During an unproblematic pregnancy** the training with ELANEE Pelvic Floor Training Aids Phase I can strengthen the pelvic floor and prevent potential subsequent problems.
- **After giving birth you can resume the training as soon as the cervix is closed (after about 6 week postpartum period).** If in doubt, ask your doctor or midwife.
- Training with ELANEE pelvic training aids can bring about the desired success in cases of **ineffective non-invasive pelvic floor training** caused by not being able to feel the pelvic floor muscles.

- **The ELANEE Pelvic Floor Training Aids are also suitable for the treatment of weaknesses of the pelvic floor resulting from a weak bladder (load/stress incontinence).** However, bladder problems have various causes. Other causes of the existing urinary incontinence should be ruled out by your doctor and/or midwife in advance.

## 5. Functionality

The pelvic floor is a flat, elastic muscle plate which is approximately palm-sized located between the pubic bone, coccyx and the buttock. It supports the internal organs in the abdominal cavity. The urethra, vagina and anus pass through the pelvic. Pregnancy and childbirth, hormonal changes/deficiency during the climacteric period or menopause, heavy physical exertion or being overweight can lead to over-stretching or weakening of the pelvic floor. The consequences can be urinary incontinence (load/ stress incontinence) or an incorrect position of the organs in the abdomen. However, the pelvic floor can be strengthened through exercise, just like any other muscle. A strong pelvic floor can also have a positive impact on sexual pleasure.

The ELANEE Pelvic Floor Training Aids Phase I active consists of a set of four different weights that enable active pelvic floor training with the biofeedback phenomenon. The biofeedback phenomenon is based on the fact that the cone slides out of the vagina through its own weight. The pelvic floor muscles react actively to this by spontaneously and autonomously contracting repeatedly in order to hold the weight. These repeated active contractions ensure muscle development - just like training muscles with weights.

Four different weights in the kit enable you to increase active exercise intensity. The training weight is increased if the muscles can hold a cone with ease. In this way a huge increase can be achieved from the lightest to the heaviest weight. An active training with pelvic floor training aids makes sense for women of all ages. It is used during pregnancy and after birth to prevent weaknesses of the pelvic floor and facilitate regeneration. In addition, it is ideal for reducing, eliminating and preventing urinary incontinence (due to load/stress incontinence) in women, particularly during the climacteric period and menopause.

## 6. Usage

1. Thoroughly **wash** the ELANEE Pelvic Floor Training Aids prior to initial usage and ensure clean conditions (e.g. wash hands).
2. **Test** to see which cone you will use to start your workout: To do this, take the **yellow (lightest) cone** and insert it slowly, completely and deeply about 3-4 cm into the vagina as you would normally insert a tampon.
  - If you can hold the yellow cone for 10 minutes while standing and walking, you can use the next heaviest cone (blue). If not, please continue with point 3!
  - If you can also hold the blue cone for 10 minutes, take the next weight, etc.
3. You complete the training with the cones by merely increasing the carrying time until the cone being used can be held **twice daily for over 10 minutes** while standing and walking, for several days in succession. Now consciously increase your feeling of your pelvic mus-

cles by actively contracting and relaxing your muscles in alternation in order to let the cone slide out and draw in again.

4. **After a few weeks** of daily training with the ELANEE Pelvic Floor Training Aids you will have learnt to feel the pelvic floor muscles without the use of cones. However, occasional practice with the help of the cones is still recommended.
5. In order to maintain the newly developed muscles and continue the training, we recommend long-term training with the **ELANEE Pelvic Floor Training Aids - Phase II.**

**Tip:** The application of lubricant is particularly helpful when training with the ELANEE Pelvic Floor Training Aids if your vagina is dry.

## 7. Safety Safety Information

-  Thoroughly clean the ELANEE Pelvic Floor Training Aids both prior to initial usage and after every application and ensure clean conditions (e.g. wash hands).
-  Please check your training aids before initial usage. Use only your cones if they are in perfect condition. If the training aid is scratched or otherwise damaged this can cause injury.
-  In cases of vaginal infections, inflammations, injury or painful complaints in the genital area and during menstruation or a problematic pregnancy, you should avoid using the Pelvic Floor Training Aids.
-  When other complaints (such as a dropping of the pelvic organs) exist, you should have other causes than a weakness of the pelvic floor ruled out in advance by a doctor

before you start with the training.

-  In cases of a known hypersensitivity and allergies to one of the named ingredients, you should not use the aids or only use them after consultation with a doctor.
-  If you are a virgin, you should not use the pelvic floor training aids.
-  The wearing of cervical caps, contraceptive rings, diaphragms or intrauterine devices during training are possible. However, you must check the correct positioning before and after each training session. If in doubt only use after consultation with a doctor.
-  Please consult your doctor if your symptoms worsen or do not improve after 6 to 8 weeks of active training.
-  Do not use steam sterilisation, dry heat sterilisation or boiling to clean the aids as this can lead to distortion of the material.

## **8. Storage and service life**

There are no special considerations.

## **9. Cleaning and disposal**

Clean the training aids after each use with warm water and soap or disinfect the cone with a commercially available skin disinfectant.

Before being used again, any possible soap residue or disinfectants should be thoroughly rinsed off the training aids with water.

The ELANEE Pelvic Floor Training Aids Phase I active can be

disposed of as standard domestic waste.

Please separate the individual components of the packaging into groups of identical materials and dispose of them in the appropriate recycling system.

When doing this please note the local and regional waste disposal regulations.

## 10. Material composition

The ELANEE Pelvic Floor Training Aids Phase I active consist of:

- an aluminium or steel core coated with polypropylene
- a polyamide yarn coated with polyurethane

## 11. Explanation of symbols used

	Manufacturer		Caution!
	Please refer to instructions		
<b>REF</b>	Article number	<b>LOT</b>	Batch number

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